

BLOOD SIMULANT CONCENTRATE

When diluted in water, TrueClot® Blood Simulant Concentrate will create a highly realistic simulated blood for use in traumatic bleeding control, wound packing and advanced hemostatic dressings training. The simulated blood will form realistic simulated blood clots when used with TrueClot® Simulated Hemostatic Gauze or TrueClot® Clotting Solution.

FEATURES:

- Small form factor, dilute with water when ready to use
- Forms highly realistic, simulated blood clots when mixed with TrueClot® Clotting Solution or Simulated Hemostatic Gauze
- Washes off of skin with soap and water (for best washability use bar soap)
- Washes out of clothing with household laundering
- 5 year shelf life
- Available in 1 gallon (3.8 L) and 4 x 1 liter sizes

DIMENSIONS:

- 1 Gallon (3.8 L) packet
Item # 5115
Dims/Weights: 6" x 9" x 1.5" (15.2 x 22.3 x 3.8 cm)
0.4 lbs (190 g)
- 4 x 1 L Packets
Item # 5117
NSN: 6910-01-727-9867
Dims/Weights: 6" x 9" x 1.5" (15.2 x 22.3 x 3.8 cm)
0.5 lbs (225 g)



Item # 5115:
1 Gallon (3.8 L) packet



Item # 5117:
4 x 1 L packets

