

## BLOOD SIMULANT PRE-MIX

TrueClot® Clotting Solution is mixed with TrueClot Blood Simulant to instantly create realistic simulated blood clots or avulsed tissue for use in moulage, bleeding control, trauma or surgical training. It may also be applied to training gauze to cause simulated clotting when packed into flowing Blood Simulant during wound packing training.

### FEATURES:

- Forms highly realistic, simulated blood clots when mixed with TrueClot® Clotting Solution or Simulated Hemostatic Gauze
- Washes off of skin with soap & water (best with bar soap)
- Washes out of clothing with household laundering
- 3 year shelf life
- Available in 1 gallon (3.8 L) and half gallon (1.9 L) volumes
- Kinematic Viscosity: 10-17 cSt at room temperature (human whole blood: 5-6 cSt)



**Item # 5112:**  
1 Gallon (3.8 L)



### DIMENSIONS:

- Blood Simulant Gallon  
Item # 5112  
NSN: 6910-01-694-8337  
Dims/Weights: 6.1" diameter x 11.6" tall  
9 lbs (15.5 x 29.5 cm - 4.1 kg)