

CLOTTING SOLUTION

TrueClot® Clotting Solution is mixed with TrueClot Blood Simulant to instantly create realistic simulated blood clots or avulsed tissue for use in moulage, bleeding control, trauma or surgical training. It may also be applied to training gauze to cause simulated clotting when packed into flowing Blood Simulant during wound packing training.

FEATURES:

- Forms highly realistic, simulated blood clots when mixed with TrueClot® Clotting Solution or Simulated Hemostatic Gauze
- Washes off of skin with soap & water (best with bar soap)
- Washes out of clothing with household laundering
- 3 year shelf life
- Powdered form (only outside of N. America) user adds distilled water
- Kinematic Viscosity: 10-17 cSt at room temperature (human whole blood: 5-6 cSt)

DIMENSIONS:

- Clotting Solution
Item # 5132
16.9 oz (500 ml)
Dims/Weights: 6" x 11.5" x 3" (15.2 x 29.2 x 7.6 cm)
1.25 lbs (567 g)
- Clotting Solution Powder
Item # 5133
Dims/Weights: 6" x 6" x 0.25" (15.2 x 15.2 x 0.6 cm)
1.3 oz (36 g) Makes 500 ml solution



Item # 5132:
Solution



Item # 5133:
Powder
(Available outside of N. America only)

