

TOURNIQUET APPLICATION TRAINER (TAT)

The TrueClot® Tourniquet Application Trainer is a wearable prop that allows full application of a tourniquet for training purposes. Simulated blood is pumped through the device to simulate hemorrhage. When a tourniquet is properly applied bleeding is stopped. The padded arm cuff reduces pinching and pain during classroom or scenario training.

FEATURES:

- Wearable, actively bleeding wound with simulated vessel (can be worn by an actor or placed on a mankin)
- Must apply a tourniquet correctly and at full force to stop bleeding
- Apply tourniquet high and tight or 2 - 3 inches above the wound
- Heavy duty padding to reduce pain
- Compatible with windlass, ratchet and other twisting-style tourniquets
- Durable ballistic nylon construction
- Easy to clean with water and air dry
- Available in Dark and Light skin tones

**Tourniquet Not Included*

DIMENSIONS:

Length: 13" (33 cm)

Circumference at Biceps: 19" (48.2 cm)

Kit in Case: 11" x 13.5" x 5.2" (28 x 34.3 x 13.2 cm)
3.2 lbs (1.5 kg)

Liquid Capacity: 33.8 oz (1 liter)



Item # 5171

NSN: 6910-01-697-8286

***NOTE:** Add "L" for Light Skin and "D" for Dark Skin tone.

